



SURGICAL,

CATALOG
2019-2020

CONTENTS

Sr #.	TITLE	Pages
01	ELEVATORS	03-010
02	ROOT TIP PICKS	11-13
03	PERIOTOMES	14
04	PERIOSTEALS	15-18
05	SURGICAL CURETTES (SPOON)	19-20
06	BONE CHISELS/ FILES / RONGEURS	21-25
07	SCISSORS	26-29
08	TUNGSTON CARBIDE TIPPED SCISSORS	30-31
09	TISSUE PLIERS / FORCEPS / NEEDLE HOLDER	32-37
10	HAEMOSTATIC FORCEPS	38-39
11	RETRACTORS	40-43
12	SCALPEL HANDLES / ASPIRATORS / SYRINGES	44-45-46
13	MISCELLANEOUS / ATRAUMATIC FORCEPS	47-48
14	EXTRACTION FORCEPS	49-59

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

MacMillan



12

E12
422111M
MacMillan
Gouge**E11M**
422312M
MacMillan
Gouge**E12M**
422413M
MacMillan
Gouge**E13M**
4225

19

E19
4228

20

E20
4229

Friedman



23

E23
422231F
Friedman
[Craig]
Gouge Serrated**E31F**
422632F
Friedman
[Craig]
Gouge Serrated**E32F**
422738
Serrated**E38**
423039
Serrated**E38**
423144C
Bradley**E44C**
42328
Crane Pick**E8**
4220

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.




46

E46
4233

46R

E46R
423469W
Woodward Hand
Pressure Chisel
E69W
423577
Serrated**E77**
423677R
Serrated**E77R**
4237

80

E80
4238

81

E81
4239

81A

E81A
4240

81E

E81E
424192
Serrated**E92**
424293
Serrated**E93**
4243

98

E98
4244

190

E190
4245

191

E191
4246

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Cogswell

A
Cogswell**EA**
4247B
Cogswell**EB**
4248

Stout

11
Stout**E11**
424911A
Stout**E11A**
4250

Heidbrink

1
Heidbrink
Straight**E1**
42512
Heidbrink
Straight**E2**
425240
Heidbrink
Straight**E40**
425341
Heidbrink
Straight**E41**
4254

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Cryer:

Cryer's right and left elevators with triangular sharp-pointed blades and hollow handles.



51
Curtis

E51
4255

52
Curtis

E52
4256



44
Cryer Small

E44
4257



45
Cryer
Small

E45
4258



25
Cryer
Medium

E25
4259



26
Cryer
Medium

E26
4260



A30
Cryer

EA30
4261



A31
Cryer

EA31
4262



31
Cryer Large

EA31
4263



32
Cryer Large

EA32
4264

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Seldin



34S
Seldin

E34S
4265



1L
Seldin

E1L
4266



1R
Seldin

E1R
4267



4L
Seldin

E4L
4268



4R
Seldin

E4R
4269



34
Seldin

E34
4270



25
Seldin

E25
4271



3S
Seldin

E3S
4272



301W
Seldin

E301W
4273



304W
Seldin

E304W
4274

Warwick James:

Small-bladed, left, right, and straight elevators.



Warwick-
James
Straight

EWJS
4275



Warwick-
James
Left

EWJL
4276



Warwick-
James
Right

EWJR
4277

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Apexo



71
Miller Apexo

E71
4278

72
Miller Apexo


E72
4279

73
Miller Apexo

E73
4280

74
Miller Apexo

E74
4281



303
Apexo

E303
4282

301
Apexo

E301
4283

4
Schmeckebier
Apexo

E4
4284

5
Schmeckebier
Apexo

E5
4285

Potts



6
Potts

E6
4286

7
Potts

E7
4287

Woodward



15
Woodward

E15
4288

16
Woodward

E16
4289

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Luxating Elevator: The luxating elevators have thin, sharp blades that are used to cut periodontal ligaments. These instruments are used by rocking and cutting verses and elevator that lifts and pries. The importance of luxating elevator is to reduce pressure on adjacent tooth.

Hourigan

H5
Hourigan
Serrated**EH5**
4290H6
Hourigan
Serrated**EH6**
4291

Luxating Elevator

Straight
2mm**EL2S**
4292Straight
3mm**EL3C**
4293Straight
5mm**EL5S**
4294Curved
2mm**EL2C**
4295Curved
3mm**EL3C**
4296Curved
5mm**EL5C**
4297


Coupland Gouge

1C
Coupland
Gouge**E1C**
42982C
Coupland
Gouge**E2C**
42993C
Coupland
Gouge**E3C**
4300A1C
Coupland
Gouge**EA1C**
4301A2C
Coupland
Gouge**EA2C**
4302A3C
Coupland
Gouge**EA3C**
4303

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

T-bar Handle

The Hourigan T-bar handle is twice the length of the Potts T-bar handle.



Cross-bar Handle




These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Heidbrink




Apical:

Set of three elevators—right, left and straight.




These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

West



Howard

Slightly heavier design than Heidbrink picks.




These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Davis Teasers



Thin sharp blades used to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament. May be used wherever an extraction is indicated, but are particularly important when considering the placement of a dental implant.

Single-Ended Periotomes



These instruments are used for retracting periosteum after removing gingival tissue.



Bennett

4400



P-149

4401



Hopkins-1

4402



Hopkins-2

4403



Mead-3

4404



Mead-3Mini

4405



Molt-9

4406

Tin Coating
Molt-9

4407



GF-14

4408

These instruments are used for retracting periosteum after removing gingival tissue.



3/9A



4409



PR-3

4410

Tin Coating
PR-3

4411

P16



4412

P15



4413

Howard-9



4414

Howard-9
Mini

4415

Allen-9

4416

These instruments are used for retracting periosteum after removing gingival tissue.



24G

4417

Tin Coating
24G

4418



Benque

4419



Buser

4420

Tin Coating
Buser

4421



KN-1

4422



KN-152

4423



Hourigan-2

4424

These instruments are used for retracting periosteum after removing gingival tissue.



HF-20

4425

Tin Coating
HF-20

4426



Molt-5

4427



Selden-23

4428



Freer

4429

OHL-10
Periosteal
Elevator


4430

Woodson-1
Modify

4431

These instruments are used for curettage and cyst removal and debridement of tooth socket.

Lucas-Curettes



These instruments are used for curettage and cyst removal and debridement of tooth socket.

Moller-Curettes

8
Miller

4485

9
Miller

4486

10
Miller

4487

11
Miller

4488

12
Miller

4489

Molt-Curettes

Molt
2/4

4490

Bone Files: For final smoothing of bone by pushing out.

Bone Chisels: For splitting and removing bone by pulling out.



1
Gardner

4540



2
Gardner

4541



5
Gardner

4542



5R
Gardner

4543



6
Gardner

4544



9
Gardner

4545



1
Chandler

4546



2
Chandler

4547



4
Chandler

4548



5
Chandler

4549



7
Chandler

4550



8
Chandler

4551



1
Wakefield

4552




2
Wakefield

4553




18
McFarland

4554




235g
Mallet

4555



130g
Mallet

4556



Mead-2
Mallet

4557

Bone Files: For final smoothing of bone by pushing out.

Miller Collburn-Cross Cut



The small paddle shaped end is half the size of the larger paddle shaped end.

1X
Miller-Collburn
Cross-Cut
4600



Curved working end combined with a straight working end.

2X
Miller-Collburn
Cross-Cut
4601



Curved working end combined with a small paddle shaped working end.

3X
Miller-Collburn
Cross-Cut
4602

Miller-Straight Cut



Same design as the 1X Cross-Cut but with Straight-Cut.

21
Miller
4603



Same design as the 2X Cross-Cut but with Straight-Cut.

45
Miller
4604

Bone Files: For final smoothing of bone by pushing out.

Miller-Straight Cut



52
Miller

4605



64
Miller

4606

Wahl-Straight Cut



1
Wahl

4607



2
Wahl

4608

Howard-Straight Cut



12
Howard

4609



67
Howard

4610

Martensite stainless steel Rongeurs are used for tooth extraction as a whole by trimming and recontouring alveolar bone. For gross tissue removal.



Blumenthal
15cm

4660



Blumenthal
45° 14.5cm

4661



Blumenthal
90° 14.5cm

4662



Mead 1A
14cm

4663



Beyer
18cm

4664



Beyer
Mini 15.5cm

4665



Luer
15cm

4666



Friedman
14cm

4667

Used for trimming tissue or cutting sutures.

Wagner



5
Wagner
Straight
11.5cm/4½"

S5

4720



6
Wagner
Curved
11.5cm/4½"

S6

4721



7
Wagner
Angled
11.5cm/4½"

S7

4722

Kelly



1
Kelly
Straight/Pointer
16cm/6¼"

S1


4723



2
Kelly
Straight/Pointer
16cm/6¼"

S2


4724



Kelly
Curved
18cm/7"

S1L

4725



Kelly
Straight
18cm/7"


S2L

4726

Martensite stainless steel Rongeurs are used for tooth extraction as a whole by trimming and recontouring alveolar bone. For gross tissue removal.




Mini-Friedman
15
14cm
4668



Forcep
Kerison
3mm
4669

Forcep
Kerison
3mm
4670



Micro Friedman
15
14.5cm
4671




Micro Friedman
90
14.5cm
4672



Micro Friedman
45
14.5cm
4673

Used for trimming tissue or cutting sutures.



Used for trimming tissue or cutting sutures.

General Surgical Scissors



21
Straight/
Pointed
14.5cm/5³/₄"

S21
4744



22
Straight/
Blunt-Sharp
14.5cm/5³/₄"

S22
4745



23
Curved/
Pointed
14.5cm/5³/₄"

S23
4746



Lister
Bandage
14.5cm/5³/₄"

SLB
4747

Suture Scissors



13
Suture
15cm/6"

S13
4748



13S
Suture
15cm/4³/₄"

S13S
4749



Crown &
Gold
Straight
15cm/4³/₄"

SCGS
4750



Crown &
Gold
Curved
15cm/4³/₄"

SCGC
4751



Crown &
Gold
Curved
Pedodontic
15cm/3.5"

SCGCP
4752

Super-Cut Scissors: The razor-sharp cutting edge cuts effortlessly and softly also through close material. The serrated cutting edge on the second blade prevents a slipping from tissue and sutures. Super-Cut is not practical for wires, cartilages and other hard material. As identification, one handle ring is gold plated.

Super cut scissors



10
Double-Curved
13.5cm/5 1/4"

S10SC
4800



14
Double-Curved
11.5cm/5 1/4"

S145C
4801



16
Goldman-Fox
12.5cm/5"

S16SC
4802



17
Iris
Straight/Delicate
11.5cm/4 1/2"

S17SC
4803



18
Iris
Curved/Delicate
11.5cm/4 1/2"

S18SC
4804



Metzenbaum
Straight/Blunt
14.5cm/5 1/2"

S5054
4805




Metzenbaum
Curved/Blunt
14.5cm/5 1/2"

S5055
4806



Metzenbaum
Straight/
Pointed
14.5cm/5 1/2"

S5056
4807



Metzenbaum
Curved/
Pointed
14.5cm/5 1/2"

S5057
4808



Metzenbaum
Delicate-
Curved/
Blunt
14.5cm/5 1/2"

S5069
4809

Used for trimming tissue or cutting sutures.





These pliers are used to grasp and stabilize soft tissue flaps during suturing.

Semkin-Taylor



31
Semkin-Taylor
Straight
12.5cm/5"

SLB
4870



32
Semkin-Taylor
Curved
12.5cm/5"

SLB
4871



33
Semkin-Taylor
1x2 Straight
12.5cm/5"

SLB
4872



34
Semkin-Taylor
1x2 Curved
12.5cm/5"

SLB
4873



16
3x4
15cm/6"

SLB
4874



50-L
Fox Locking 3x3
15cm/6"


SLB
4875



Kramer-Nevins
Tissue Pliers
1x2 Curved

SLB
4876

Suture Pliers



20 Corn
Suture Pliers
15cm/6"

SLB
4877



41
Adson Plain
12cm/4 3/4"

SLB
4878



42
Adson 1x2
12cm/4 3/4"

SLB
4879



43
Adson-Brown
7x7
12cm/4 3/4"

SLB
4880



Russian
15cm/6"

SLB
4881



Russian
20cm/8"

SLB
4882



Gerald Plain
18cm/7"
Straight/Serrated

SLB
4883



Gerald 1x2
15cm/7"
Straight

SLB
4884

Allison



Allison
Baby 4x5
12.5cm/5"

SLB
4885



Allison
Straight 5x6
15cm/6"

SLB
4886

Adson

Adson
12cm/4 $\frac{3}{4}$ "**TCTPA**
4887Adson 1x2
12cm/4 $\frac{3}{4}$ "**TCTPA**
4888Adson
Brown
12cm/4 $\frac{3}{4}$ "**TCTPAB**
4889

General

General
15cm/6"**SLB**
4890General
1x2
15cm/6"**SLB**
4891

Gerald

Gerald
18cm/7"**SLB**
4892Gerald 1x2
18cm/7"**SLB**
4893Gushing
18cm/7"**SLB**
4894

Cushing

Cushing
Angled
18cm/7"**SLB**
4895Cushing
Bayone
19cm/7 $\frac{1}{2}$ "**SLB**
4896Potts Smith
18cm/7"
(Shown)**SLB**
4897Lister
Bandage
14.5cm/5 $\frac{3}{4}$ "**SLB**
4898De Bakey
15cm/6"
(Shown)**SLB**
4899De Bakey
19cm/7 $\frac{1}{2}$ "**SLB**
4900

Martensite stainless steel instruments are used for grasp and guide the needle while suturing.




Hegar-Baumgartner
(Gardner 5)
14cm/5 1/2"

MHHB
4950



Webster
Smooth
13cm/5"

NHW
4951



Mayo-Hegar
(Gardner 6)
Grooved
16cm/6 1/4"

NHMH
4952



Crile-Wood
Grooved
15cm/6"

NHCW
4953



Olsen-Hegar
Needle
Holder/Scissors
17cm/6 3/4"

NHOH
4954



Boynton
(Hu-Friedy)
12cm/4 1/4"

NHB
4955



Castroviejo
14cm/5 1/2"

NHCV
4956



Collier
12.5cm/5"

NHC
4957




Derf
Grooved
12.5cm/5"

NHD
4958



Mathieu
(Machatt)
14cm/5 1/2"

NHM
4959



Arruga
Serrated
16cm/6 1/2"

NH5022
4970



Webster
Smooth
11.5cm/4 1/2"

NH5030
4971



Crile-Wood
15cm/6"

NH5038
4972



Mayo-Hegar
16cm/6 1/4"

NH5042
4976

Olsen-Hegar (Needle Holder)



Olsen-Hegar
Needle Holder
Scissors
17cm/6 3/4"

NH5068
4973



Olsen-Hegar
Needle Holder
Scissors
14cm/5 1/2"

NH5072
4974



DeBakey
18cm/7"
Delicate

NH5046
4975



Mayo-Hegar
18cm/7"

NH5044
4977

Used for grasp and guide the needle while suturing.



Matheu-Kocher
14cm/5½"

NH5074
4978



Matheu
14cm/5½"

NH5076
4979



HMT
Thin Jaw
14cm/5½"

NH5050
4980



Fine Swedish
Pattern
15cm/6"

NH5052
4981



Micro-Vascular
Pattern
15cm/6"

NH5082
4982



Micro-Vascular
Pattern
18cm/7"

NH5084
4983



Mini-Ryder
Pattern
15cm/6"

NH5092
4984



Mini-Ryder
Pattern
18cm/7"

NH5094
4985

Used for grasp and guide the needle when suturing.



Derf
11.5cm/4½"

NH5032
4960

Hegar-
Baumgartner
11.5cm/4½"

NH5034
4961

Halsey



Halsey
Smooth
13cm/5"

NH5036
4962

Halsey
13cm/5"

NH5037
4963

Castroviejo



Castroviejo
Straight
14cm/5½"

NH5020
4964

Castroviejo
Curved
14cm/5½"

NH5021
4965

Micro
Castroviejo
Straight
14cm/5½"

NH5020M
4966

Micro
Castroviejo
Curved
14cm/5½"

NH5021M
4967

Castroviejo
Straight
18cm/7"

NH5024
4968

Castroviejo
Straight
Round Handle
18cm/7"

NH5024R
4969

A multipurpose instrument used to clamp off blood vessels, remove small root tips and grasp loose objects. All Hemostats

Hartman



Hartman
Straight
10cm/4"

HHS
5030



Hartman
Curved
10cm/4"

HHC
5031



5
1x2 Teeth
Straight
12.5cm/5"

HHS
5034

Adson



Adson
Straight
19cm/7"

HAS
5032



Adson
Curved
19cm/7"

HAC
5033



Rochester-pean
Curved
16cm/6"

HRPC
5035

A multipurpose instrument used to clamp off blood vessels, remove small root tips and grasp loose objects. All Hemostats

Halsted-Mosquito



3
Halsted-Mosquito
Curved
12cm/4"

HHS
5036



3
Halsted-Mosquito
Straight
12cm/4"

HHS
5037



Kocher
1x2 Teeth Curved
14cm/5 1/2"

HKC
5038

Kelly



1 Kelly
14cm/5"

HHS
5039



2 Kelly
Curved
14cm/5"

HHS
5040



Kelly-Rankin
Straight
16cm/6"

HHS
5041



Kelly-Rankin
Curved
16cm/6 1/4"

HHS
5042



Skin Hooks: For hooking skin during surgery.

Retractors: For holding back the lips and cheeks.




Skin Hooks: For hooking skin during surgery.

Retractors: For holding back the lips and cheeks.





Retractors: For holding back the lips and cheeks, and tongue.

Mouth Gag: To prop open the mouth and prevent it from closing.



Retractors: For holding back the lips and cheeks, and tongue.



Scalpel Handle: For holding a mess blade.

	3 12cm
5200	
	3G 12cm
5201	
	4 13.5cm
5202	
	5 14.5cm
5203	
	5A 14.5cm
5204	
	Double Head 14.5cm
5205	
	7 16.5cm
5206	
	3L 21cm
5207	
	Scalpel 15.2cm
5208	
	Scalpel 180° 16.2cm
5209	

These tools are used for precision section at surgical sites.
Tools shown are 50% size.

Shown 50% size.



Actual tip size (Internal diameter)

5250	NH5052	○	1 Coupland Tip	1.5mm
5251	NH5052	○	2 Coupland Tip	2.0mm
5252	NH5052	○	3 Coupland Tip	2.5mm
5253	NH5052	○	4 Coupland Tip	3.0mm

Actual tip size (Internal diameter)

5254	NH5052	○	1 Coupland Tip/Side Holes	1.5mm
5255	NH5052	○	2 Coupland Tip/Side Holes	2.0mm
5256	NH5052	○	3 Coupland Tip/Side Holes	2.5mm
5257	NH5052	○	4 Coupland Tip/Side Holes	3.0mm

Shown 50% size.



Byrd
Self-Cleaning
20cm/8"

NH5074
5258

Actual tip size (Internal diameter)

5259	NH5052	○	B1 Byrd	1.5mm
5260	NH5052	○	B2 Byrd	2.0mm
5261	NH5052	○	B3 Byrd	2.5mm


Shown 50% size.



Byrd
Self-Cleaning
20cm/8"

NH5074
5262

Shown 50% size.



Actual tip size (Internal diameter)

5263	NH5052	○	6 Frazier	1.5mm
5264	NH5052	○	6 Frazier	2.0mm
5265	NH5052	○	10 Frazier	2.5mm
5266	NH5052	○	12 Frazier	3.0mm

Frazier
Suction Tube

Oral Surgery
22cm/8.5"
1mm

NH5074
5267

Yandauer
26.5cm/10.5"

NH5074
5268

Andrew-Pychon
23cm/9"

NH5074
5269

Shown 50% size.

Shown 50% size.


Shown 50% size.




Locking Pliers
16.5cm/6 1/2"

NH5074
5270

Syringes: An instrument for injecting aesthetics for a cartridge.
3-Way Syringes - Supplies air, water, or a mist through separate interval channels created by combining pressurized air with water flow.




Intraligamentary
Syringe
1.8ml
5320



Syringe
1.8ml
5322




Syringe
1.8ml
5324



Syringe
1.8ml
5326




Paroject
5328




Syringe
1.8ml
5321



Syringe
1.8ml
5323



Syringe
1.8ml
5325



Syringe
2.2ml
5327



Paroject
5329

Miscellaneous tools for laboratory.



1
Lwanson
Spring Caliper
CLP1
5380



2
Lwanson
Spring Caliper
CLP2
5381



Bone
Morselizer
BM
5382



Tofflemire
Junior
TJ
5383



Tofflemire
Senior
TS
5384




Tofflemire
Universal
TU
5385



Boley
Gauge
BG
5386

Our Unique tip design makes it mold to the individual shape of your patient's tooth. The smooth handle with openings, allows for more thorough and easier cleaning. The newly styled handle allows a more uniform transmission while applying pressure, pulling and for rotary movements. More comfortable design can be used by right or left-handed practitioners! Our handles are lighter and are easier to use.



Upper
Universal

5430




Lower
Universal

5431




Upper
Anterior

5432




Lower
Anterior

5433




Upper
3rd Molar

5434




Lower Molar
Modify
Cow-Horn
Bead

5435




Lower Root
Extra Long

5436



Upper Root

5437



Lower
Wisdoms

5438



Lower
Roots

5439




Upper Roots

5440

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.


Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




5
Lower Third
Molars

5441




15
Lower
Molars

5442




16S
Lower Molars

5443



17
Lower Molars

5444



16
Lower Molars
Cow-Horn Beak

5445



1AF
Upper
Anteriors

5446




Upper Roots

5447




Upper Roots

5448



Splitting
Forces

5449




Upper Centrals
Laterals &
Cuspids

5450

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.


Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




99C
Upper
Premolars
& Incisors

5451




88L
Upper Molars
Left Side

5452




88L
Upper Molars
Right Side

5453




10S
Upper Molars
Straight-Handle

5454




53L
Upper Molars
Left Side

5455




32
Upper Molars
& Premolars
Strong Beak

5456




23
Lower Molars
Cow-Horn Beak

5457




18R
Upper Molars
Right Side

5458




69
Fragments
& Small Roots
Upper & Lower

5459



65I
Upper Incisors
& Roots

5460



Universal
for Upper
Teeth & Roots

5461




53R
Upper Molars
Right Side

5462

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.


Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




151A
Lower
Premolars
Incisors & Roots

5463




151XAS Children
Lower Premolars
Incisors &
Roots Canines

5464




151AS
Split Lower
Premolars
Incisors &
Roots Canines

5465




151AF Apical
Beak Lower
Premolars
Incisors &
Roots Canines

5466




101
Universal
for Premolars

5467




150
Upper Premolars
Incisors &
Roots Canines

5468




150AF
Apical Beak
Upper Premolars
Incisors &
Roots Canines

5469




150AS
Split Beak
Upper Premolars
Incisors &
Roots Canines

5470




150XAS
Children Split
Upper Premolars
Incisors &
Roots Canines

5471



150A
Upper Premolars
Incisors &
Roots Canines

5472



150S
Upper Premolars
Incisors &
Roots for
Children

5473




151
Upper Premolars
Incisors &
Roots

5474

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.


Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




151S
Lower Premolars
Incisors &
Roots for Children

5475




210H
Upper Wisdoms
with Hook Handle

5476




210S Upper
Wisdoms
with Straight
Handle

5477




222
Lower Wisdoms

5478



222AF
Apical Beak
Lower Wisdoms

5479




286
Upper Roots
Premolars
& Incisors

5480

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




67 AF
Apical Beak
Upper Wisdoms

5493




68
Lower Roots

5494



73
Lower Molars
(Hawks Bill)

5495



74
Lower Incisors
Roots

5496



31
Lower Roots

5497



32
Lower Molars

5498




33
Lower Roots

5499




33M
Lower Roots
With Long Beak

5500




67
Upper Wisdoms

5501




51A
Upper Roots
Narrow Beak

5502



45
Lower Roots

5503



44
Upper Roots

5504

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



5481



5482



5483



5484



5485



5486



5487



5488



5489



5490



5491



5492

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




95
Upper Molars
Left

5505




121
Upper Wisdoms

5506



167 Universal
Patter for
Lower Teeth
Molars & Premolars

5507




168
Upper
Canines

5508




74N
Lower Incisors
& Roots
Canines

5509



74NF
Lower Incisors
& Roots

5510



79
Lower Wisdoms

5511



79A
Lower Wisdoms

5512



86C
Lower Molars

5513



89
Upper Molars
Right

5514



90
Upper Molars
Left

5515



94
Upper Molars
Right

5516

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




F-1
Upper
Universal

5517




F-2
Lower
Universal

5518




F-3
Upper Anterior

5519




F-4N
Lower Anterior

5520




169
Lower Cnines

5521




233
Lower Roots
Parallel Beaks

5522



235
Upper Roots
Parallel Beaks

5523




MD-1
Universal
Upper Wisdom
Premolars &
Roots

5524




MD-4
Lower Molars

5525



22L
Lower Wisdoms
Left

5526




22R
Lower Wisdoms
Right

5527

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.


SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




F-7C
Lower Molar
Modify
Cow-Horn Beak
5528




F-8
Upper Universal
5529




F-9
Lower Universal
5530



F-10
Lower Molar
5531



F-11
Upper Molar
Left
5532




F-12
Upper Molar
Right
5533




F-13
Upper Molar
Right
5534



F-14
Upper Molar
Left
5535




F-4
Lower Anterior
5536



F-5
Lower Anterior
5537



F-6
Upper
3rd Molars
5538



F-7
Lower
3rd Molars
5539



F-1
Upper Universal
Pedo

5551



F-2
Lower Universal
Pedo

5552



F-3
Upper Anterior
Pedo

5553



F-4
Lower Anterior
Pedo

5554



F-5
Lower Anterior
Pedo

5555



F-6
Upper 3rd
Pedo

5556



F-6
Lower 3rd
Pedo

5557



F-16
Lower Root
Pedo

5558



F-17
Lower Root
Pedo

5559



B23
Lower Molar
Pedo

5560



150-SK
Upper Universal

5561



151-SK
Lower Universal

5562




Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.


SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.




F-15
Upper
Premolars
Universal
5540



F-16
Upper
Root
5541




F-17
Lower Root
5542




F-18
Upper Root
Fragment Narrow
5543




F-18TN
Upper Roots
Fragment Narrow
5544




F-19
Lower Roots
Fragment Narrow
5545




F-19NW
Lower Roots
Fragment Narrow
5546



46XL
Lower Root
Extra Long
5547



F-20
Upper
Premolar
1ST & 2RD Molars
5548



49
Upper Roots
5549



51XL
Upper Roots
5550



+92 524 360 300



+92 333 332 2214



www.aaazinstruments.com



info@aaazinstruments.com



1 KM Head Marala Road
Murad Pur, Sialkot, Pakistan.